



KIDDRC  
Histology Services  
Protocols

For more info contact: Jing Huang, Histology Specialist  
913-588-5996

## Basic Immunocytochemistry procedure (lab method)

- 1.** Deparaffinize slides and run down through graded alcohol to dHOH.  
For frozen sections: Let slides come up to room temperature and air dry.  
If fixed, hydrate slides in 3 changes of dHOH.  
If unfixed, fix in cold acetone for 10 minutes.  
Let slides air dry. Rehydrate with three changes of dHOH. If using another fixative, fix slides for 10-30 minutes, and wash in three changes of dHOH.
- 2.** 3% HOOH for 15 minutes.
- 3.** Rinse twice with DHOH, then once with TBS.
- 4.** Incubate slides with 10% Normal Goat Serum (NGS) for 20 minutes.
- 5.** Incubate slides with primary antibody for appropriate time.  
overnight at 4°C  
2hrs. at room temp.  
1 hr. at room temp. or 37°C
- 6.** Rinse 3 with TBS (Rinse X3 quick, then 2-5 minutes washes.)
- 7.** Incubate with biotinylated GAR or GAM 1:200 for 30 min. at room temp.
- 8.** Rinse X 3 with TBS.
- 9.** Incubate with HRP-avidin D (or streptavidin D) 1:1000 for 10 min. at room temp.
- 10.** Rinse X2 with TBS. Rinse X1 with distilled water (not deionized/deionized water can inhibit DAB)

**11.** Develop slides---approx. 5min. at room temp. Observe reaction under microscope.

DAB-Zymed Kit: 1 ml distilled water

1 drop buffer concentrate rgnt. (bottle 1)---mix

1 drop rgtn 2 (DAB chromogen 20 X)--mix

1 drop rgnt 3 (HOOH)---mix

After developing, rinse in 3 changes of distilled water, then use DAB enhancer, if desired. Add 1 drop of rgnt. 4 to 1 ml of distilled water and mix.

Cover tissue and incubate 1-3 minutes at room temp. Color development may be monitored under scope.

Rinse in 2-3 changes of water.

DAB-Vector Kit: 2.5 ml. distilled water

1 drop buffer sol'n--mix

2 drops DAB, mix

1 drop HOOH, mix

(can add 1 drop nickel sol'n for black chromogen)

After development, rinse, then incubate slides in copland jar with 1% Copper Sulfate for 5 minutes. Rinse and procede.

True Blue (Kirkegard-Perry): Let come to Room temp, use full strength.

\* Primary Ab should be at a higher dilution (lower concentration) than with DAB.

**12.** Dehydrate slides to xylene and coverslip with a nonaqueous mountant.

### **To Use Alkaline Phosphatase:**

**1.** After secondary Ab, use Alkaline Phosphatase-avidin D as enzyme 1:50 for 20 min. at room temp.

**2.** Rinse X3 with TBS.

**3.** Develop with Vector Red: usually about 10-20 minutes at room temp.

2.5 ml 100 mM Tris HCL pH 8.2

1 drop Rgt. 1, mix

1 drop Rgt.2, mix

1 drip Rgt.3, mix

1 drop levamisole

\*Vector Red is also fluorescent with rhodamine, fluorescein or AMCA excitation filter.

**4.** Dehydrate to xylene and coverslip.

### **TBS(TRIS BUFFER SALINE---0.5M)**

Stock solution .5M

9.875g Tris Base

65.94g Tris HCL

fill to 1 liter with distilled water

80g NaCl

pH to 7.6

Working solution

100ml stock TBS

900ml distilled water

### **Dilution Buffer**

8.640 ml TBS

960  $\mu$ l Normal goat serum

400  $\mu$ l 1% BSA

10  $\mu$ l .02% sodium azide

use this buffer for diluting the primary antibody

Sodium azide: for 25 ml of .02%

0.005g/25ml

100mM Tris HCl (for alkaline phosphatase)

1.32g Tris Base

2.215g Tris HCl

250 ml d water

pH to 8.2