

Equipment



IITC Cold Analgesiometer

Used in assays of:

Cutaneous Sensitivity to Cooling

Cold Allodynia

Thermal (Cold) Analgesia

One IITC incremental cold plate is available for use with rats or mice. The unit is capable of latency or threshold based sensory testing. The plate is microprocessor-controlled, and can be used for incremental warming or cooling from any starting temperature from 0 deg. C to 70 deg. Centigrade in 0.1 degree increments. The unit can also double as a constant temperature warm or cold analgesiometer.