

Equipment



Accuscan Thermal Analgesia Monitor

Used in assays of:

**Cutaneous Sensitivity to Thermal Stimulation
Thermal Hyperalgesia
Thermal Analgesia**

The RBF has an Analgesia Monitor (model # AI0407AM0207) consisting of a thermally controlled plate, which is used to measure complex sensory responsiveness to warming. The controller holds the plate at the set temperature indefinitely without over-shooting. The unit also includes an electronic timer, which accurately measures the response latency of the subject. General features include a sophisticated temperature controller (digital set point, auto tuner), built-in electronic timer with 1 ms resolution, a time palm switch for ease of operation, and an acrylic arena cover.